


# Development of an Electro-Oculography (EOG) Measurement System

Christopher Fabian, Mike Fuller, Bin Guo,  
Xiefan Lin and Mike Kavanagh

BIOM 611 Project Presentation  
April 30, 2002

## Electro-oculogram

- Measurement of corneal-retinal potential
- Three electrode setup for horizontal eye movement detection




## EOG Overview of Function

- Eye position:
  - Left, potential is positive
  - Right, potential is negative
  - Straight, EOG output is zero
- Approximately linear over 30 degree arc
- Accurate to 1 degree

## Applications

- Sleep and dream research
- Reading ability and visual fatigue




La sorcière transforma  
la petite souris en princesse.  
Elle désirait lui faire épouser  
le prince merveilleux.

© Metrovision

## Applications – Cont'd

- Retinal dysfunction
  - Some disorders of the retina exhibit abnormal or absent changes in the corneal-retinal potential during dark and light adaptation processes.
- Vestibular and balance dysfunction
  - Nystagmus: characteristic slow-phase/fast-phase eye movement measurable by EOG
  - prevent drunk drivers from taking to the road.

## Circuit Design of EOG




3-op-amp differential amplifier

3-pole Butterworth lowpass filter

Inverting amplifier

## Actual EOG Circuit




## Testing


- Left, right, left, right
- Left, center, right, center
- Visual pursuit
  - Sinusoid @ 0.1 Hz
  - Sinusoid @ 0.4 Hz
- Optokinetics
  - Rotating light bars
  - Evokes nystagmus



## Left, Center, Right, Center




## Left, Right, Left, Right



## Visual Pursuit @ 0.1 Hz




## Visual Pursuit @ 0.4 Hz



## Brief Description of Nystagmus


Nystagmus is an involuntary eye movement characterized by a slow-phase and fast-phase.

### Dix-Hallpike Right



## Nystagmus – Cont'd

- ENG can be used to diagnose balance and dizziness problems
- Possible causes of nystagmus:
  - Vestibular-ocular reflex
  - Nervous system dysfunction (e.g. alcohol)
  - Visual Field stimulus
    - Optokinetics



© Micromedical Technologies

# Optokinetics



# Optokinetics – Cont'd





## Test Result Summary

- Overall, worked very well!
- Measurement system exhibited linear behavior within the  $\pm 30^\circ$  range of motion
- Drift was an issue, but resolved with periodic calibration
- Noise reduced while preserving required bandwidth

## Cleaning up the Signal

- Calibration
  - The corneo-retinal potential is not stable and varies with a number of factors, in particular the level of light. Light source should be adjustable by the use of filters or other means.
- Anti-drift
  - DSP, calibration procedure, etc.

## Cleaning up the Signal – Cont'd

- A masking cone
  - May be used to reduce electrical interference (actually provided by some manufacturers).
- Artifact rejection
  - Eliminates some of the obvious peaks or drifts in the distorted recorded waveforms due to blinks and other movement.

## Making It Clinically Usable

- This simple design provides virtually good results, so there is not much left to do with the circuit itself.
- A/D conversion & Computer Program:
  - Artifact rejection
  - Robust control of lights.
    - Usually alternately illuminated for 1 second each, ten times per minute.
  - Waveform display & user interface
  - Data analysis
  - Decision-making procedures

## DISCLAIMER:

No group members were (irreparably)  
harmed during the “testing” of this “device.”

## Acknowledgment

We would like to thank Jorge Gonzalez, M.S.  
of the  
UVA Vestibular and Balance Center  
for his assistance in the testing of our instrument.

Questions?